

* The casing design is subject to change.



Handheld Vision Screener
easyRef Lite

probably the fastest objective vision screener


- Objective and accurate refractive assessment
- Ideal for challenging patients such as infants, disabled patient
- Revolutionary design brings superior user experience

EASY OPERATION THAT YOU CAN MASTER IN SECONDS



Step 1

Simply put the screener against the forehead



Step 2

Centre the cross on the pupil



Step 3

Save, print or share the report

Sphere range	-10.0D ~ +10.0D
Sphere tolerance	±0.5D
Cylinder range	-5.0D ~ +5.0D
Cylinder tolerance	±0.5D
Axis of astigmatism range	1° ~ 180°
Data acquisition time	1s/eye
External	External must connect to PC and IPAD
Net weight	360g

Specification subject to change without notice

MOPTIM[®]

Technology Research and Development:
Shenzhen Moptim Imaging Technique Co.,Ltd.
www.moptim.com
sales@moptim.cn

Manufacturer:
Shenzhen Certainn Technology Co.,Ltd.
www.certainn.com
sales@moptim.cn

Address:
Bldg. 2-C, Section 2, GOTO Digital Technology Park,
No.137 Bulan Rd., Buji Subdistrict, Longgang District,
Shenzhen 518112, P.R. China